

IHE Integration Statement			
Vendor	Product Name	Version	Date
Canon Medical Informatics	Vitrea Read	9.1	2024-12-03

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Access to Radiology Information	Image Display	None
Access to Radiology Information	Report Reader	None
Consistent Presentation of Images	Evidence Creator	None
Consistent Presentation of Images	Image Display	None
Chest X-Ray CAD Display	Image Display	None
Digital Breast Tomosynthesis	Image Display	None
Key Image Note	Evidence Creator	None
Key Image Note	Image Display	None
Mammography Image	Image Display	None
Portable Data for Imaging	Image Display	None
Portable Data for Imaging	Portable Media Importer	None
Portable Data for Imaging	Report Reader	None
Simple Image and Numeric Report	Report Reader	None

Internet address for vendor's IHE information :

<https://www.mi.medical.canon/support-and-services/customer-support/ethics-and-compliance/compliance-documents/>

Links to Standards Conformance Statements for the Implementation

HL7	https://www.mi.medical.canon/support-and-services/customer-support/ethics-and-compliance/compliance-documents/
DICOM	https://www.mi.medical.canon/support-and-services/customer-support/ethics-and-compliance/compliance-documents/

Link to general information on IHE : www.ihe.net